

SUGGESTED SAMPLE QUANTITIES

Homogeneous powdered products (i.e., starch, flour, flavors)	50 g
Coarse heterogeneous material (i.e., snack mix, chips)	100 g
Whole seeds or kernels/beans	2000 g
Other miscellaneous products (i.e., leaves)	usually 50 - 100 g
Lecithin	50 g
Liquid Products	50ml

- All samples should be packed securely in unbreakable material (plastic or paper materials) to avoid contamination during shipping.
- Samples sent in unstructured packaging (i.e., plastic or paper bagging) should be sealed in a second plastic bag.
- Liquid samples should be sent in leak-proof containers with the lids/tops sealed and which are then placed in a plastic bag in case of leaking or breakage.